

CERTIFICATE OF ANALYSIS

MONOCLONAL ANTIBODIES

Human Thyroid Stimulating Hormone (TSH) Mouse Monoclonal Antibody

Code No.	BM216-V6D5
Grade	Affinity Purified
Lot No.	

RECEIVER INFORMATION

Expiry Date	
Storage Temperature	2-8°C
Storage Notes	Avoid repeat freeze thaw cycles
Shipping Temperature	2-8°C

PRODUCT INFORMATION

Specificity	Antibody specifically detects Human TSH beta subunit. Minimal cross-reactivity with TSH alpha subunit, LH, FSH or hCG.
BBI Solutions Clone No.	V6D5
Application	Lateral Flow. Can be used in Lateral Flow as the Detection antibody pairing with BM216-S1A5 as the Capture antibody.
Isotype	IgG ₁
Host	Mouse
Purification	Protein A
Presentation	Single homogenous batch, supplied liquid in 0.9% NaCl and 0.1% Sodium Azide.

HEALTH AND SAFETY

Application	For in-vitro use and research or further processing only.
Precaution	The product has been labelled T (toxic) under UK regulations: COSHH Sch 1(6) for sodium azide concentration $\geq 0.1\%$.

ANALYSIS

TESTS	SPECIFICATIONS	RESULTS
-------	----------------	---------

Determination Method:

Protein Concentration by
Optical Density at 280nm using
 $E_{0.1\%,1\text{cm}} = 1.4$

≥ 0.5 mg/ml

mg/ml

Name:		Position:	
Signed:		Date:	