

CERTIFICATE OF ANALYSIS

MONOCLONAL ANTIBODIES

NS1 Dengue virus, subtype 1, mouse monoclonal antibody

Code no.	BM409-G3D1
Grade	Affinity Purified
Lot no.	892-201

RECEIVER INFORMATION

Expiry Date	2 nd August 2016
Storage Temperature	2-8 °C (up to 1 month) below -15°C (long term)
Storage Notes	Avoid repeat freeze thaw cycles
Shipping Temperature	2-8°C

PRODUCT INFORMATION

Specificity	Antibody specifically detects NS1, human Dengue virus. (Subtype 1).
Immunogen	Purified from Native Dengue virus, NS1.
Source	Ascites of balb/c mouse.
BBI Solutions clone no.	G3D1
Application	Recommended for use in immunoassays: ELISA, Western Blot and Lateral Flow assays.
Isotype	IgG
Host	Mouse
Purification	Protein G
Presentation	Single homogenous batch, supplied liquid in a PBS buffer, pH7.2 containing 0.1% sodium azide.
Purity	≥ 90%.

HEALTH AND SAFETY

Application	For in-vitro use and research or further processing only.
Precaution	The product has been labelled T (toxic) under UK regulations: COSHH Sch 1(6) for sodium azide concentration ≥ 0.1%.

ANALYSIS

TESTS	SPECIFICATIONS	RESULTS
-------	----------------	---------

Determination method:

Protein Concentration by
Optical Density at 280nm using
 $E_{1\%,1\text{cm}} = 1.38$

>0.5 mg/ml

1.38 mg/ml

Purity

Defined by a major band $\geq 90\%$ on coomassie stained SDS-PAGE.

Name:		Position:	
--------------	--	------------------	--

Signed:	Date:
----------------	--------------